

Zygomatic²

PLUS

SCIENTIFIC JOURNAL

IMPLANTAT

Discover IMPLANTAT,
the educational habitat of S.I.N.
An online teaching platform created to make
more professionals accelerate their career
and increase their success.

Access
IMPLANTAT.GLOBAL
or scan the QR CODE
and begin your journey of
knowledge now!



EDUCATION POWERED BY S.I.N.



Zygomatic²

PLUS

SCIENTIFIC JOURNAL

AUTHORS:

DR. FÁBIO JOSÉ BARBOSA BEZERRA

Specialist in Implantology by INEPO – SP; Master in Periodontics by UNIP – SP; PhD student in Biotechnology at IBB – UNESP Botucatu. Author of 5 books, 18 book chapters and over 60 articles published on scientific journals. Bezerra has also lectured in 16 countries until now.

DRA. BRUNA GHIRALDINI

PhD. in Implantology by the Paulista University, MBA in Business Management – Getúlio Vargas Foundation. Master in Periodontics – Paulista University, Specialist in Periodontics – Paulista University. Post-Graduate in Implantology – INEPO. Post-graduate in Aesthetics – SOESP. Coordinator of the Product Research and Development – S.I.N.

DR. FERNANDO GIOVANELLA

Graduated in Dentistry by the University of Vale do Itajaí; Specialized in Implantology by FURB-SC and Oral Surgery by ABP-PR; Specialized in Oral and Maxillofacial Surgery and Traumatology by UFPel-RS; Master in Implantology by the Holy Heart University – Bauru/SP; Professor of Anesthesiology and Oral and Maxillofacial Surgery at the Regional University of Blumenau (FURB).

DR. BRENNNO MARCONDES NEGRI

Specialist in Implantology by UNIP – SP; Specialist in Public Health and Family Dentistry by Unyleya – SP; Graduated at FO-USP; Consulting Scientist at the Product Research and Development Department – S.I.N..

DR. GUSTAVO PETRILLI

Graduated in dentistry by the University of São Paulo; Specialist in Implantology by INEPO – São Paulo; Specialist in Oral-Maxillofacial Traumatology and Surgery by UNIBAN – São Paulo; Scientific Consultant at S.I.N.

DR. FABIO AZEVEDO

Graduated in Dentistry by the University of Santo Amaro (UNISA); Specialist in Implantology by SENAC; Master in periodontics by UNG; MBA in Health and Clinic Management by FGV; Specialist consultant at the R&D department at S.I.N.

TABLE OF CONTENTS


05 • EDITORIAL

06 • HISTORY & DEVELOPMENT

08 • CLINICAL CASES

10 New Zygomatic Implants for Resolving Conventional Implant Failures and Maxillary Sinus Lift, by Dr. Fernando Giovanella and Dr. Mauricio Motta

17 New Zygomatic Implants for the Resolution of Atrophic Maxilla, by Dr. Christopher Barrett and Dr. Rajan Sheth



EDITORIAL

Zygomatic implants have gained scientific reputation through studies that prove the applied technique's success, both in the short and long term. Currently, specialists consider this solution an alternative in the therapeutic portfolio, making it possible to achieve the best results in different challenging and complex cases of atrophic maxilla.

ZYGOMATIC EXPLORES THE MAXIMUM EFFICIENCY AND PRECISION IN THE MOST CHALLENGING SITUATIONS IN DENTAL PRACTICE!

Despite its popularity, the installation of zygomatic implants demands a high level of knowledge since there are still dangers and potential complications during the procedure.

With the advent of new digital tools, an evolution in surgical techniques has been seen. Implants started to be positioned with the head, in many cases, on top of the crest of the alveolar ridge. On some occasions, by looking at the work model or clinically in the mouth, it is impossible to differentiate between a well-positioned zygomatic implant and a well-installed conventional implant.

Given that, planned and assertive perforations have drastically reduced the risk of mistakes and deviations in perforations that are potentially harmful to noble anatomical structures.

Sinus problems have been reduced with the possibility of exteriorizing the zygomatic implant, with the proper positioning of the implant's head, always supported by a residual bone bed.

The evolution of the surgical process has been clear and beneficial to patients; for that reason, the development of new geometries for zygomatic implants has been necessary throughout the years.

By: DR. FERNANDO GIOVANELLA

HISTORY AND DEVELOPMENT

DEVELOPMENT

The initial design with coils and a wider head, ideally developed for the Branemark Technique, is still frequently used. Recently, some changes have been adopted for improvement, meeting the clinical demands within the contemporary notion of correct indications for the proper use of the technique.

CLINICAL INDICATIONS

The zygomatic implant is recommended for complex and challenging situations that arise in the implantologist's practice. It is indicated for severe cases of atrophic maxilla, where the zygomatic implant head is often supported by a bone in the rim, with the implant's body and head sometimes without bone coverage.

At this point, we reach the challenge of the current technique: peri-implant tissue stability in the long term.

PRECAUTIONS

It is not unusual to find well-installed and rehabilitated zygomatic implants; however, they may exhibit tissue dehiscence and exposure of some cervical coils. This dehiscence can cause concern for the patient, making hygiene harder to perform and potentially leading to inflammations or even local infections. Therefore, the search for an optimized implant design for the clinical situation is essential for the solution and surgical procedure of the cases.

Considering this, the contemporary design of a zygomatic implant must take into account osseointegration and primary stability, mainly in the implant's apical area. Additionally, it must present a tissue-friendly design in the cervical area. This means it should have the smallest volume possible at the cervical area and feature a macrostructure that does not irritate soft tissues. Regarding the surface, a specific treatment is something that will further optimize stability and tissue adherence.

HANANO SURFACE

The inclusion of a bioactive surface elevates the product to the next level, increasing predictability and clinical success. The HANano surface is composed of hydroxyapatite nanocrystals, whose size and form are similar to those of human bone. It has a 20-nanometer width, sintered onto titanium, which promotes a change in the surface energy. In addition, it increases hydrophilicity and provides a substrate that stimulates higher protein adsorption, adherence, proliferation, and differentiation of osteoblasts.

Preliminary studies have shown excellent fibroblast adherence, which is highly desirable in the cervical area. Considering these factors, the new Zygomatic by S.I.N. intends to be a solution in line with modern and updated indications for this type of implant, as well as enabling higher stability, easier osseointegration, and promoting good peri-implant health.

CONCLUSION

Although the External Hexagon is the most studied prosthetic platform regarding zygomatic implants, the use of internal connections has gained the interest of many dental surgeons. The fact that these do not demand the use of assemblers, along with the possibility of personalizing angulations with mini-pillars, are perceived as advantages by professionals. Therefore, having two prosthetic platform options in the portfolio makes sense to serve the global market.


In summary, the technique of zygomatic implants is an excellent treatment alternative, with vast scientific evidence. Its surgical technique has evolved, as has the implant's design, which has followed the development of the technique.

The new Zygomatic by S.I.N. brings characteristics in accordance with the clinical challenge, combined with good practical usability. Its diverse components for prosthetic solutions and proven treatment provide the implantologist with the best zygomatic implant on the market.



Fernando Giovanella, PhD.



CLINICAL CASES



NEW ZYGOMATIC IMPLANTS TO SOLVE THE UNSUCCESS OF CONVENTIONAL IMPLANTS AND MAXILLARY SINUS LIFT



GUSTAVO PETRILLI, MASTER, & FÁBIO AZEVEDO, MASTER.

Gustavo Petrilli, Graduated in dentistry by the University of São Paulo; Specialist in Implantology by INEPO – São Paulo; Specialist in Oral-Maxillofacial Traumatology and Surgery by UNIBAN – São Paulo; Scientific Consultant at S.I.N.

Fábio Azevedo, Graduated in Dentistry by the University of Santo Amaro (UNISA); Specialist in Implantology by SENAC; Master in periodontics by UNG; MBA in Health and Clinic Management by FGV; Specialist consultant at the R&D department at S.I.N.

CASE INTRODUCTION

The case begins with a series of failed attempts at surgeries with conventional implants and maxillary sinus lifts. Given this unsuccessful trajectory, the proposed approach to fully rehabilitate the arch involved the implementation of four zygomatic implants. This sequence of events highlights the complexity of the clinical condition and the need for an innovative solution for the patient’s oral restoration.

TREATMENT EXPLANATION


The patient had a history of failed conventional implant surgeries. Subsequently, she underwent a maxillary sinus lift and the placement of new implants, but this treatment was also unsuccessful. Given this situation, the use of four zygomatic implants was proposed and implemented to fully rehabilitate the arch.

PATIENT FILE

- Gender:** Female.
- Age:** 70 years old.
- Complaint:** Masticatory inability, lack of confidence in social relations, unstable total prosthesis, and frustration with previous treatments.
- Anamnesis:** Normotensive patient, with no underlying pathologies, in good health.
- Planning:** Installation of four zygomatic implants to rehabilitate the entire upper arch.


BEFORE

Patient with a remaining implant, severely atrophic maxilla, with history of two treatments with conventional implants and graft, both unsuccessful.




RADIOGRAPHY/TOMOGRAPHY

The first radiography shows the presence of 4 installed implants, but insufficient for a total rehabilitation. At the second radiography, it is possible to see the addition of 4 more implants in a procedure of maxillary sinus lift, which was unfortunately unsuccessful.




RADIOGRAPHY/TOMOGRAPHY

The tomography show the unsuccessful grafting and the loss of installed implants.




RADIOGRAPHY/TOMOGRAPHY

Tomography of the zygomatic bone used to plan and conduct the case.




STEP BY STEP PROCEDURE

In the intraoral view and during total detachment, it is possible to observe a severe atrophy in height and width, indicating the need for zygomatic fixations with precision.




STEP BY STEP PROCEDURE

Full-window approach, with careful detachment of the sinus membrane and access to the zygomatic body.




INTRAOPERATIVE

Installation of implants on the right side with close-up view, keeping the integrity of the sinus membrane.




INTRAOPERATIVE

Occlusal view with positioned Mini Abutments, presenting excellent emergency and prosthetic positioning.




INTRAOPERATIVE

Final photograph of the surgery with performed suture, preserving soft tissues and keeping a healthy gum phenotype.



IMPLANT USED

Close-up view of the zygomatic implant that was installed.




PROSTHETIC PHASE

Bimaxillary protocol prosthesis installed, with intraoral photograph of the case finished in 48 hours.



FINAL RADIOGRAPHY

Case finished in immediate load in less than 48 hours, presenting passivity and full settling of the structure.



NEW ZYGOMATIC IMPLANTS TO SOLVE THE UNSUCCESS OF CONVENTIONAL IMPLANTS AND MAXILLARY SINUS LIFT



GUSTAVO PETRILLI, MASTER, & FÁBIO AZEVEDO, MASTER.

Gustavo Petrilli, Graduated in dentistry by the University of São Paulo; Specialist in Implantology by INEPO – São Paulo; Specialist in Oral-Maxillofacial Traumatology and Surgery by UNIBAN – São Paulo; Scientific Consultant at S.I.N.

Fábio Azevedo, Graduated in Dentistry by the University of Santo Amaro (UNISA); Specialist in Implantology by SENAC; Master in periodontics by UNG; MBA in Health and Clinic Management by FGV; Specialist consultant at the R&D department at S.I.N.

CASE INTRODUCTION

The patient, who had previously undergone an unsuccessful implant procedure, chose to undergo a new procedure involving the use of four zygomatic implants. This approach was adopted with the intention of fully rehabilitating the dental arch.

TREATMENT EXPLANATION


The patient had a history of failed conventional implant surgeries. Consequently, the use of four zygomatic implants was proposed and implemented to fully rehabilitate the arch.

PATIENT FILE

Gender: Female.
Age: 50 years old.
Complaint: Unstable total prosthesis and frustration with previous treatments. The patient desired fixed teeth and improved social life.
Anamnesis: Normotensive patient, with no underlying conditions, in good health.
Planning: Installation of four zygomatic implants to rehabilitate the entire upper arch.


BEFORE

Patient with remaining unsatisfactory implants, severely atrophic maxilla and history of treatment with conventional implants, both unsuccessful.




RADIOGRAPHY/TOMOGRAPHY

The images show the presence of two installed implants that remain but are insufficient for total rehabilitation.




RADIOGRAPHY/TOMOGRAPHY

Tomographic maxilla cuts show the lack of bone availability for the use of conventional implants. It is an extremely atrophic case.




RADIOGRAPHY/TOMOGRAPHY

Tomography of the zygomatic bone used to plan and conduct the case.




STEP BY STEP PROCEDURE

Initial intraoral view, along with the total detachment conducted, reveals a severe atrophy in height and width. The surviving implants were removed, indicating the need to apply the zygomatic fixation technique.




INTRAOPERATIVE

The implants installed on the right and left side, respectively, are in close contact with the bone bed, with settling over the rim.




INTRAOPERATIVE

Occlusal view with mini abutments positioned, presenting excellent emergency and prosthetic positioning.




INTRAOPERATIVE

Final photograph of the surgery with the suture performed, preserving soft tissues and keeping a good gum phenotype.



IMPLANT USED

Close-up view of the Zygomatic implant that was installed.



PROSTHETIC PHASE

Bimaxillary protocol prosthesis installed. Intraoral photograph of the case finished in 48 hours.



FINAL RADIOGRAPHY

Case finished in immediate load in less than 48 hours, with full passivity and settling of the structure.



OUR GLOBAL PRESENCE



POINT YOUR PHONE'S CAMERA AT THE
QR CODE AND SEE WHERE S.I.N. IS PRESENT

HEADQUARTERS

2140 Vereador Abel Ferreira Av
Jardim Anália Franco
São Paulo – SP - Brazil

FACTORY

421 Soldado Ocimar Guimarães
da Silva St - Jardim Anália Franco
São Paulo – SP - Brazil

PORTUGAL BRANCH

General Ferreira Martins St, 10 8D -
1495-137 Algés - Portugal
+351 21 412-0336

sinimplantsystem.com

SUPERIOR QUALITY AND TECHNOLOGY

*WE WARRANT, BECAUSE WE ARE
PROUD OF OUR PRODUCTS.*



S.I.N.'s main priority is assuring the quality and safety to our clients.
Offering the best for implants, components, surgical kits and tooling is the
base of all our action.

INSPECTION IN A 100% OF THE BATCHES MANUFACTURED

The quality control is made in all S.I.N. products, to assure the success in the
surgeries of all our clients, to meet the best quality standards, as well as to
add value to all the ones who chose to give a smile back to people.



**IMPLANTS WITH WARRANTY
FOR LIFE***



**5 YEARS OF WARRANTY
PROSTHESIS COMPONENTS***



*SCAN THE LATERAL QR CODE
TO ACCESS S.I.N WARRANTY
TERMS OR ACCESS THE LINK
sinimplantsystem.com/lifetime-warranty/





Get to know Implantat, the educational streaming of S.I.N.

 [implantat.global](https://www.implantat.global)

internacional@sinimplantsystem.com
www.sinimplantsystem.com

Visit our social networks



[@sinimplantglobal](https://www.instagram.com/sinimplantglobal)



[@sinimplantglobal](https://www.facebook.com/sinimplantglobal)



[S.I.N. Implant System](https://www.youtube.com/S.I.N.ImplantSystem)



[/sin_implant](https://www.tiktok.com/sin_implant)