


M A N U A L

# BEYOND Full Arch



# IMPLANTAT

EDUCATION POWERED BY S.I.N.

Discover IMPLANTAT, the educational habitat of S.I.N.  
An online teaching platform created to make more professionals  
accelerate their career and increase their success.




Access  
**IMPLANTAT.GLOBAL**  
or scan the QR Code  
and begin your journey  
of knowledge now!



|                                                     |           |
|-----------------------------------------------------|-----------|
| <b>BEYOND FULL ARCH.....</b>                        | <b>4</b>  |
| <b>IDENTIFYING THE PATIENT.....</b>                 | <b>6</b>  |
| <b>COMPLETE PLANNING.....</b>                       | <b>8</b>  |
| <b>TREATMENT GUIDELINES.....</b>                    | <b>9</b>  |
| SURGICAL PLANNING.....                              | 9         |
| AXIAL POSITIONING - 4 TO 8 IMPLANTS.....            | 9         |
| ANGLED POSITIONING - 4 TO 8 IMPLANTS.....           | 11        |
| <b>S.I.N. PRODUCTS.....</b>                         | <b>14</b> |
| <b>EPIKUT S.....</b>                                | <b>15</b> |
| EPIKUT S IMPLANTS.....                              | 15        |
| INSTALLATION .....                                  | 16        |
| KIT.....                                            | 16        |
| MILLING.....                                        | 17        |
| <b>EPIKUT S LONG.....</b>                           | <b>19</b> |
| EPIKUT S LONG IMPLANTS.....                         | 19        |
| INSTALLATION.....                                   | 20        |
| KIT.....                                            | 21        |
| PREPARATION .....                                   | 22        |
| MILLING.....                                        | 22        |
| <b>ZYGOMATIC.....</b>                               | <b>24</b> |
| ZYGOMATIC IMPLANTS.....                             | 24        |
| INSTALLATION.....                                   | 26        |
| KIT.....                                            | 28        |
| MILLING.....                                        | 29        |
| <b>S.I.N. SOLUTIONS.....</b>                        | <b>32</b> |
| DIGITAL IMPRESSIONS .....                           | 33        |
| INTRAORAL SCANNER.....                              | 34        |
| PROSTHETIC REHABILITATION.....                      | 35        |
| <b>CONSIDERATIONS ON THE FULL ARCH CONCEPT.....</b> | <b>36</b> |
| <b>SCIENTIFIC PUBLICATIONS .....</b>                | <b>37</b> |

# BEYOND Full Arch



**DARE TO GO BEYOND**

Dare to go beyond is to offer your patients complete protocols that restore essential daily functions, such as eating and smiling, while boosting self-esteem and transforming lives, even in severe cases of maxillary atrophy.

### **Dare to go beyond, explore the art of Beyond Full Arch!**

The total loss of natural teeth is typically seen in the elderly, requiring a rehabilitation that fully restores function and aesthetics. To meet these expectations precisely and reliably, **S.I.N.** presents **Beyond Full Arch**, which is based on a Decision Tree protocol. This protocol is suitable for patients with different degrees of bone resorption in the maxilla and mandible, whether mild, moderate, or severe.

The goal of this treatment is to restore both the smile and chewing function of patients.

The concept of full arch rehabilitation involves the use of various types of implants, varying in diameters, lengths, and techniques, to ensure the best personalized solution for each case.

In this manual, you will find all the **Beyond Full Arch** solutions that **S.I.N.** offers for different clinical indications.

# IDENTIFYING THE PATIENT

Some preoperative exams may be necessary, including:

- › Radiographic exam
- › Cone Beam Computed Tomography (CBCT)
- › 3D Software
- › Surgical guides
- › Reverse planning

Some measures are necessary to better identify a patient qualified for a full rehabilitation procedure, as well as to improve the treatment hierarchy:

› **01** Any medical conditions that could affect the patient's treatment outcome or their suitability for surgery should be considered. The patient's expectations should also be analyzed and documented.

› **02** The patient's expectations and their history of implant failures should be analyzed, as well as any parafunctional habits, such as teeth grinding and bruxism.

› **03** The initial radiographic analysis with a panoramic X-ray is important. However, the professional should perform a computed tomography or a cone beam tomography. For full arch rehabilitation, this exam is mandatory for proper planning.

› **04** Evaluate the condition of the remaining teeth. In many cases, some teeth may need to be removed for prosthetic reasons. The evaluation should consider the state of periodontal health, untreatable focal infections, as well as the condition of the patient's soft tissues and mucosa. In the context of treatment longevity, mucosal conditions will be important.

› **05** The 3D planning software can be an essential and significant tool for better planning and achieving a better final result.

› **06** When clinically possible, the use of surgical guides and 3D planning using specialized software can be an important ally in improving results and accuracy.

› **07** Reverse planning: Ideally, all planning and surgical treatment should follow a proper prosthetic evaluation and pre-planning.

# COMPLETE PLANNING

The complete arch planning should be entirely focused on the prosthesis. From this base, we begin an analysis of the quantity and quality of bone.

- 1** - Residual hard and soft tissues.
- 2** - Transition line: to determine aesthetic and prosthetic options.
- 3** - Evaluation of the patient's residual bone and degree of bone resorption: to determine surgical and prosthetic procedures.
- 4** - Digital planning:

3D models based on the patient's radiological data, along with radiographic guidance, allow the professional to assess the quality and quantity of available bone. Vital anatomical structures, such as the alveolar nerve and maxillary sinus, can also be marked, so that prosthetic planning can be carried out.


For a successful full rehabilitation, the initial prosthetic planning, built with the help of a tomographic guide and according to the correct occlusion using Cone Beam Computed Tomography (CBCT) images, will define the proper implant position, considering anatomical structures, especially for tilted posterior implants. If the patient already has a well-established conventional full denture, it can also be used as a guide and as an immediate provisional prosthesis supported by implants. Additionally, specialized planning software can be used to determine the position and angulation of the implants.

## BONE RESORPTION PATTERNS

**SEVERE  
RESORPTION**



**PARTIAL  
RESORPTION**



**NORMAL  
VOLUME**



**THE RESIDUAL ALVEOLAR BONE DIRECTS THE SURGICAL PROTOCOL.**

# TREATMENT GUIDELINES

Ideally, a full arch rehabilitation on implants should optimize immediate loading and implant installation with a minimum torque of 40N.cm. Once a full arch treatment plan has been developed and approved, the procedure and the full arch concept can be implemented.

## SURGICAL PLANNING AXIAL POSITIONING - 4 TO 8 IMPLANTS


In cases of skeletal class III, eight implants should be considered for installation in an axial or tilted position. In most cases, four to six implants are ideal.



Full arch with four implants.



Full arch with six implants.



Full arch with eight implants.

- Clin Oral Implants Res. 1995 Dec;6(4):227-31.

**Ten-year survival rates of fixed prostheses on four or six implants ad modum Brånemark in full edentulism.**  
Brånemark PI, Svensson B, van Steenberghe D.

- Implant Dent. 2015 Dec;24(6):680-5.

**Stress Distribution in Bone and Implants in Mandibular 6-Implant-Supported Cantilevered Fixed Prosthesis: A 3D Finite Element Study.**


Padhye OV, Herekar M, Patil V, Mulani S, Sethi M, Fernandes A.

## UPPER ARCH


The distribution of implants should consider the shape and size of the maxillary sinus and regions with less available bone.

Make an incision that allows full access to the alveolar ridge, followed by a full flap. In case of tooth extractions, careful curettage of the socket is required, as well as the creation of a bone plateau for implant placement. The transition of the smile line and lip support should be considered.


Implants should be installed axially in the anatomical positions of the central incisor, the canine, the first premolar, and the first molar. Anterior implants should be placed palatally, and the installation of implants at 45° angle posteriorly should be considered, depending on the patient's anatomy. Distal cantilevers beyond one tooth and severe parafunction should be avoided. For immediate loading, a primary implant stability of >45Ncm is recommended. In areas of tooth extraction, implants should preferably be placed in the interalveolar regions.



Installation of straight anterior implants (maxilla).



Installation of posterior implants (maxilla).



Installation of straight abutments (maxilla).



Instalação de abutments retos (maxila).

## LOWER ARCH

The distribution of implants should consider the position and distance of the inferior alveolar and mental nerves.

Make an incision allowing full access to the alveolar ridge, followed by a full flap. In the case of tooth extractions, careful curettage of the socket is required, as well as the creation of a bone plateau for implant placement. The transition of the smile line and lip support should be considered, as well as the location of the mental foramen, which will be the distal limit of the implants, especially if there is insufficient bone height in the posterior regions of the mandible.


Implants should be installed axially in the anatomical positions of the central incisor, the canine, the first premolar, and the first molar. Anterior implants should be placed lingually, and the installation of posterior implants with an inclination of 45 degrees should be considered, depending on the patient's anatomy. Distal cantilevers beyond one tooth and severe parafunction should be avoided. For immediate loading, a primary implant stability of >45N.cm is recommended. In extraction areas, implants should preferably be placed in the interalveolar regions.



Installation of straight anterior implants (mandible).



Installation of posterior implants.



Installation of straight abutments (mandible).




Full arch, fully installed (mandible).

## ANGLED POSITIONING - 4 TO 8 IMPLANTS



Four implants (tilted distal implants)

To treat cases with a high degree of bone resorption, it is possible to use four implants: two anterior axial implants and two posterior implants tilted up to 45°.



Six implants (tilted distal implants)

The inclination of the two posterior implants improves the bone-implant contact and provides excellent bone support, even with minimal bone volume. A surgical guide, considering the patient's anatomy, can be used if necessary.

- Clin Implant Dent Relat Res 2018 Oct;20(5):867-874  
**Implant rehabilitation of the edentulous jaws: Does tilting of posterior implants at an angle greater than 45 affect bone resorption and implant success?: A retrospective study.**  
Malchiodi L, Moro T, Cattina DP, Cucchi A, Ghensi P, Nocini PF
- J Oral Maxillofac Surg 1999; 57(3):281-287  
**Implant treatment without bone grafting in severely resorbed edentulous maxillae.**  
Mattsson T, Köndell PA, Gynther GW, Fredholm U, Bolin A
- Int J Oral Maxillofac Implants 2009 May-Jun;24(3):527-33  
**Palatal Positioning of Implants in Severely Resorbed Edentulous Maxillae.**  
Peñarrocha M, Carrillo C, Boronat A, Balaguer J, Peñarrocha M
- Clin Implant Dent Relat Res. 2005;7 Suppl 1:S88-94  
**All-on-4 immediate-function concept with Brånemark System implants for completely edentulous maxillae: a 1-year retrospective clinical study.**  
Maló P., Rangert B., Nobre M.

- Clin Implant Dent Relat Res. 2014 Dec;16(6):836-55  
**The all-on-four treatment concept: a systematic review.**  
Patzelt SB, Bahat O, Reynolds MA, Strub JR.
- J Prosthodont Res. 2017 Apr;61(2):123-132.  
**Biomechanical analysis of immediately loaded implants according to the "All-on-Four" concept.**  
Horita S, Sugiura T, Yamamoto K, Murakami K, Imai Y, Kirita T.
- Int J Oral Maxillofac Implants 2016;31:1017-1022. (Portuguese Edition)  
**Nova Proposta Para Reabilitação de Maxila Atrófica: Implante Inclinado LONG.**  
Marcelo de Carvalho, Liliane Pacheco de Carvalho, Rogerio de Lima Romeiro, Carlos Eduardo Francischone, Bruno Salles Sotto-Maior, Fabio Bezerra
- Compend Contin Educ Dent. 2016;37(7):458-465  
**All-on-4® Implant Treatment: Common Pitfalls and Methods to Overcome Them.**

## UPPER ARCH

The distribution of implants should consider the shape and size of the maxillary sinus and regions with less bone availability.

Make an incision allowing full access to the alveolar crest, followed by a full flap. In case of tooth extractions, careful curettage of the socket is required, as well as the creation of a bone plateau for implant placement. The transition of the smile line and lip support should be considered.

Two to four anterior implants should be installed in the palatal position, and two posterior implants with an inclination of up to 45°.

The inclination of the two posterior implants improves bone-implant contact and provides excellent bone support even with low bone volume, as it allows the use of longer implants, as well as bicorticalization. Additionally, the inclination of the implants in the maxilla provides better anchorage to the anterior bone of higher quality.


Bicortical anchorage to the cortical bone of the sinus wall and the nasal fossa can also be considered. Distal cantilevers beyond one tooth and severe parafunctions should be avoided.

For immediate loading, a primary implant stability of  $>45\text{Ncm}$  is recommended. In extraction areas, implants should be placed between the extraction sockets.

### PROSTHESIS (6 IMPLANTS)



Installation of angled implants



Installation of angled abutments




Implants and prosthesis.

### PROSTHESIS (4 IMPLANTS)



Installation of angled implants



Installation of angled abutments



Implants and prosthesis.

## LOWER ARCH

The distribution of implants should consider the position and extension of the inferior alveolar and mental nerves.


Make an incision allowing full access to the alveolar ridge, followed by a full-thickness flap. In case of dental extractions, careful curettage of the socket is necessary, as well as the creation of an osseous plateau for implant insertion. Lip support should be considered.

For fixed prostheses, two axial implants should be placed near the midline, along with two distal implants with an inclination of no more than  $45^\circ$ . This allows for screw exit and positioning of the prosthetic abutment closer to the first lower molar region, reducing the prosthetic lever arm.


For immediate loading, an implant primary stability of  $>45\text{Ncm}$  is recommended. In extraction sites, implants should be placed between the extraction sockets.



Implants positioned for full rehabilitation (inclined distals)



Abutments.



Implants and prosthesis.



# PRODUCTS S.I.N.

The S.I.N. Implant System offers various implant lines that make the daily practice of professionals easier for surgical procedures. Through Beyond Full Arch, we have the Epikut S line, with conventional implant sizes, and the Epikut S Long line, which includes longer sizes, reaching up to 24 mm. There is also the Zygomatic Plus line, which features a 16° connection. Below, we will discuss the details of each line.



# EPIKUT S

## THE IMPLANTS

The Epikut S implants are made of commercially pure, biocompatible titanium (Grade IV cold worked) and are available with only the surface treatment of dual acid etching (DAA) or with the DAA treatment followed by the application of nanometric thickness hydroxyapatite (HANano).

The exclusive macrogeometry featuring progressive cutting threads makes the Epikut the state-of-the-art choice for immediate load cases, low-density bone, and post-extraction sockets. Extremely versatile, the Epikut S can also be used in other clinical situations as long as the recommended clinical drilling protocol is followed. Its finer apex, with a diameter of 2 mm, facilitates bicorticalization.

The Epikut S line has a 16° Cone Morse connection (MT 16°), with lengths of 8.5, 10, 11.5, 13, and 15 mm, and diameters of 3.5, 3.8, 4.0, 4.5, and 5.0, as described in the table below:

| TECHNICAL MEASUREMENTS |           |                                |                              |                  |
|------------------------|-----------|--------------------------------|------------------------------|------------------|
| DAA CODE               | PLUS CODE | A<br>PLATFORM DIAMETER<br>(mm) | B<br>APICAL DIAMETER<br>(mm) | C<br>LENGTH (mm) |
| ILM 3585               | ILM 3585N | 3,5                            | 2,0                          | 8,5              |
| ILM 3510               | ILM 3510N |                                |                              | 10,0             |
| ILM 3511               | ILM 3511N |                                |                              | 11,5             |
| ILM 3513               | ILM 3513N |                                |                              | 13,0             |
| ILM 3515               | ILM 3515N |                                |                              | 15,0             |
| ILM 3885               | ILM 3885N | 3,8                            | 2,0                          | 8,5              |
| ILM 3810               | ILM 3810N |                                |                              | 10,0             |
| ILM 3811               | ILM 3811N |                                |                              | 11,5             |
| ILM 3813               | ILM 3813N |                                |                              | 13,0             |
| ILM 3815               | ILM 3815N |                                |                              | 15,0             |
| ILM 4085               | ILM 4085N | 4,0                            | 2,0                          | 8,5              |
| ILM 4010               | ILM 4010N |                                |                              | 10,0             |
| ILM 4011               | ILM 4011N |                                |                              | 11,5             |
| ILM 4013               | ILM 4013N |                                |                              | 13,0             |
| ILM 4015               | ILM 4015N |                                |                              | 15,0             |
| ILM 4585               | ILM 4585N | 4,5                            | 2,95                         | 8,5              |
| ILM 4510               | ILM 4510N |                                |                              | 10,0             |
| ILM 4511               | ILM 4511N |                                |                              | 11,5             |
| ILM 4513               | ILM 4513N |                                |                              | 13,0             |
| ILM 4515               | ILM 4515N |                                |                              | 15,0             |
| ILM 5085               | ILM 5085N | 5,0                            | 2,95                         | 8,5              |
| ILM 5010               | ILM 5010N |                                |                              | 10,0             |
| ILM 5011               | ILM 5011N |                                |                              | 11,5             |
| ILM 5013               | ILM 5013N |                                |                              | 13,0             |
| ILM 5015               | ILM 5015N |                                |                              | 15,0             |




## INSTALLATION

The placement of conventional size implants from the Epikut S line follows a technical approach that prioritizes stability and bone integration. These implants are designed with dimensions that promote anchorage in native bone, reducing the need for grafts and facilitating installation in regions of high-quality bone.

## KIT

For the installation of Epikut implants, it is necessary to use the specific kit for this line. A complete and compact kit with a linear and intuitive sequence.



| CODE     | DESCRIPTION                           |
|----------|---------------------------------------|
| FL 20    | LAUNCH DRILL Ø2,0MM                   |
| FHE 27   | CONICAL DRILL Ø2,7MM                  |
| FHE 30   | CONICAL DRILL Ø 3,0MM                 |
| FHE 33   | CONICAL DRILL Ø 3,3MM                 |
| FHI 36   | CONICAL DRILL Ø 3,6MM                 |
| FHI 38   | CONICAL DRILL Ø 3,8MM                 |
| FHI 40   | CONICAL DRILL Ø 4,0MM                 |
| FHI 43   | CONICAL DRILL Ø 4,3MM                 |
| FHI 48   | CONICAL DRILL Ø 4,8MM                 |
| ID 2720L | DIRECTION INDICATOR Ø2,7xØ2,0MM LONG  |
| ID 3020L | DIRECTION INDICATOR Ø3,0xØ2,0MM LONG  |
| ID 3320L | DIRECTION INDICATOR Ø3,3xØ2,0MM LONG  |
| ID 3620L | DIRECTION INDICATOR Ø3,6xØ2,0MM LONG  |
| ID 3820L | DIRECTION INDICATOR Ø3,8x Ø2,0MM LONG |
| ID 4020L | DIRECTION INDICATOR Ø4,0xØ2,0MM LONG  |
| ID 4320L | DIRECTION INDICATOR Ø4,3xØ2,0MM LONG  |
| ID 4820L | DIRECTION INDICATOR Ø4,8xØ2,0MM LONG  |

| CODE     | DESCRIPTION                           |
|----------|---------------------------------------|
| FH 20    | HELICOIDAL DRILL Ø2,0MM               |
| FC 41    | COUNTERSINK DRILL Ø4,1MM              |
| EXFN     | EXTENSION KEY FOR DRILL               |
| CDH 1224 | DIGITAL HEX 1.2 MEDIUM KEY            |
| CBD 01   | BI-DIGITAL KEY                        |
| MTCM 11  | TRANSMUCOSAL GAUGE CM 11,5°           |
| MTCM 16  | TRANSMUCOSAL GAUGE MT 16°             |
| CTMD 24  | CONTRA-ANGLE KEY FOR IMPLANT CM LONG  |
| CTMD 20  | CONTRA-ANGLE KEY FOR IMPLANT CM SHORT |
| CCM 24   | RATCHET KEY FOR IMPLANT CM LONG       |
| CCM 20   | RATCHET KEY FOR IMPLANT CM SHORT      |
| CTWD 24  | CONTRA-ANGLE KEY FOR IMPLANT HE LONG  |
| CTWD 20  | CONTRA-ANGLE KEY FOR IMPLANT HE SHORT |
| CCW 24   | RATCHET KEY FOR IMPLANT HE LONG       |
| CCW 20   | RATCHET KEY FOR IMPLANT HE SHORT      |
| SOP 20   | DEPTH PROBE                           |
| TMECC 02 | SURGICAL TORQUE WRENCH                |

## MILLING

### FOR BONE SOFT TYPE

Drilling sequence  
used for type IV  
bone.



Epikut S Epikut S Plus

|         |                    | 1.200 RPM    | 800 RPM       |               |               |               |                |               |               |               |
|---------|--------------------|--------------|---------------|---------------|---------------|---------------|----------------|---------------|---------------|---------------|
|         | Ø<br>DIÂM.<br>(mm) | FL 20<br>(A) | FHE 27<br>(B) | FHE 30<br>(C) | FHE 33<br>(D) | FHI 36<br>(E) | FHI 38<br>(E+) | FHI 40<br>(F) | FHI 43<br>(G) | FHI 48<br>(H) |
| ILM35xx | 3,5                | •            | •             |               |               |               |                |               |               |               |
| ILM38xx | 3,8                | •            | •             | •             |               |               |                |               |               |               |
| ILM40xx | 4,0                | •            | •             | •             | •             |               |                |               |               |               |
| ILM45xx | 4,5                | •            | •             | •             | •             | •             |                |               |               |               |
| ILM50xx | 5,0                | •            | •             | •             | •             | •             | •              | •             |               |               |



# EPIKUT S LONG

## THE IMPLANTS

The Epikut S Long implants are made of commercially pure, biocompatible titanium (Grade IV cold worked) and have a surface treatment of dual acid etching (DAA) followed by the application of nanometric layer hydroxyapatite (HANano).

They are exclusively indicated for multiple rehabilitations in fully edentulous maxillae. These implants can be combined with conventional, pterygoid, and/or zygomatic implants.

According to the complete arch rehabilitation technique, at least four implants should be distributed in the maxilla to retain the full prosthesis.

The Epikut S Long line features a 16° Cone Morse connection (MT 16°), with lengths of 18, 20, 22, and 24 mm, and diameters of 3.8, 4.0, and 4.5, as described in the table below:

| TECHNICAL MEASUREMENTS |           |                        |                      |             |
|------------------------|-----------|------------------------|----------------------|-------------|
|                        |           | A                      | B                    | C           |
| CODE DAA               | CODE PLUS | PLATFORM DIAMETER (mm) | APICAL DIAMETER (mm) | LENGTH (mm) |
| ILM 3818               | ILM 3818N | 3,8                    | 2,0                  | 18          |
| ILM 3820               | ILM 3820N |                        |                      | 20          |
| ILM 3822               | ILM 3822N |                        |                      | 22          |
| ILM 3824               | ILM 3824N |                        |                      | 24          |
| ILM 4018               | ILM 4018N | 4,0                    | 2,0                  | 18          |
| ILM 4020               | ILM 4020N |                        |                      | 20          |
| ILM 4022               | ILM 4022N |                        |                      | 22          |
| ILM 4024               | ILM 4024N |                        |                      | 24          |
| ILM 4518               | ILM 4518N | 4,5                    | 2,95                 | 18          |
| ILM 4520               | ILM 4520N |                        |                      | 20          |
| ILM 4522               | ILM 4522N |                        |                      | 22          |
| ILM 4524               | ILM 4524N |                        |                      | 24          |



## INSTALLATION

The technical approach for placing long implants is similar to that of conventional-sized implants. The Epikut S Long implants have dimensions and lengths that provide excellent bone anchorage. Additionally, placing the implants in the patient's native bone helps avoid bone grafts and facilitates implant placement in regions of excellent bone quality.

The tangential approach to the maxillary sinus, with the placement of inclined implants in the canine pillar region, utilizing the entire available bone height, allows for high torques and, consequently, the use of the immediate load technique.


Four implants are used: two posterior and two anterior. The two posterior implants are placed in the position of the premolars, with an inclined direction up to a maximum of 45°, in a mesial direction tangent to the anterior wall of the maxillary sinus, with the prosthetic platform emerging in the position of the second premolar or first molar. The two anterior implants can be placed straight or inclined.

The maximum torque for installing the Epikut S Long, without causing damage or fractures, is 80 N/cm.

Thus, the Epikut S Long line allows for the rehabilitation of edentulous patients in a simpler, faster, and more predictable manner.

## KIT

For the installation of Epikut S Long implants, it is necessary to use the specific kit for this line. It is a complete and compact kit with a linear and intuitive sequence. It features angled direction indicators divided by colors, a more modern system that facilitates its clinical use.



| CODE     | DESCRIPTION                           |
|----------|---------------------------------------|
| FL 2024  | HELICOIDAL DRILL Ø 2,0 X 24MM         |
| FHE 2324 | CONICAL DRILL Ø 2,35 X 24MM           |
| FHE 3024 | CONICAL DRILL Ø 3,0 X 24MM            |
| FHI 3324 | CONICAL DRILL Ø 3,3 X 24MM            |
| FHI 3624 | CONICAL DRILL Ø 3,6 X 24MM            |
| FHI 3824 | CONICAL DRILL Ø 3,8 X 24MM            |
| FHI 4024 | CONICAL DRILL Ø 4,0 X 24MM            |
| FHI 4324 | CONICAL DRILL Ø 4,3 X 24MM            |
| CTMD 20  | CONTRA-ANGLE KEY FOR IMPLANT CM SHORT |
| CTMD 24  | CONTRA-ANGLE KEY FOR IMPLANT CM LONG  |
| CCM 20   | RATCHET KEY FOR IMPLANT CM SHORT      |
| CCM 24   | RATCHET KEY FOR IMPLANT CM LONG       |
| IDA 17   | DIRECTION INDICATOR ANGULADO 17°      |
| IDA 30   | DIRECTION INDICATOR ANGULADO 30°      |
| IDA 45   | DIRECTION INDICATOR ANGULADO 45°      |
| CBD 01   | BI-DIGITAL KEY                        |
| CDH 1220 | DIGITAL HEX 1.2 SHORT                 |
| TMECC 02 | SURGICAL TORQUE WRENCH                |
| SOPL 01  | DEPTH PROBE LONG                      |

## PREPARATION

After the installation of Epikut S Long implants, proceed with the selection of abutments, provisional prosthetic rehabilitation, and fabrication of the final prosthesis. This step is no different from the technique used with conventional implants.

## MILLING

The use of long implants is indicated to achieve good anchorage in areas with better bone density, such as the canine pillar area and the lateral wall of the nasal cavity. The planning is based on the principles that guide all rehabilitations: clinical study of radiographic and tomographic exams, laboratory tests, digital photographs, obtaining plaster models, mounting on a semi-adjustable articulator, and creating a surgical guide through reverse planning. This planning can also be done digitally, so that the implant installation can be simulated using software and prototypes.

### FOR BONE SOFT TYPE

Drilling sequence  
used for type IV  
bone.



|         |                    | 1.200 RPM      | 800 RPM         |                 |                 |                 |                  |                 |                 |
|---------|--------------------|----------------|-----------------|-----------------|-----------------|-----------------|------------------|-----------------|-----------------|
|         | Ø<br>DIAM.<br>(mm) | FL 2024<br>(A) | FHE 2324<br>(B) | FHE 3024<br>(C) | FHI 3324<br>(D) | FHI 3624<br>(E) | FHI 3824<br>(E+) | FHI 4024<br>(F) | FHI 4324<br>(G) |
| ILM38xx | 3,8                | •              | •               | •               |                 |                 |                  |                 |                 |
| ILM40xx | 4,0                | •              | •               | •               | •               |                 |                  |                 |                 |
| ILM45xx | 4,5                | •              | •               | •               | •               | •               |                  |                 |                 |

## FOR BONE MEDIUM TYPE

Drilling sequence  
used for type II and  
III bone.



Epikut S  
Longo



Epikut S  
Longo Plus

|         |                    | 1.200 RPM      | 800 RPM         |                 |                 |                 |                  |                 |                 |
|---------|--------------------|----------------|-----------------|-----------------|-----------------|-----------------|------------------|-----------------|-----------------|
|         | Ø<br>DIÂM.<br>(mm) | FL 2024<br>(A) | FHE 2324<br>(B) | FHE 3024<br>(C) | FHI 3324<br>(D) | FHI 3624<br>(E) | FHI 3824<br>(E+) | FHI 4024<br>(F) | FHI 4324<br>(G) |
| ILM38xx | 3,8                | •              | •               | •               | •               | •               | •                | •               | •               |
| ILM40xx | 4,0                | •              | •               | •               | •               | •               | •                | •               | •               |
| ILM45xx | 4,5                | •              | •               | •               | •               | •               | •                | •               | •               |

• Use of the drill is optional.

## FOR BONE HARD TYPE

Drilling sequence  
used for type I bone.



Epikut S  
Longo



Epikut S  
Longo Plus

|         |                    | 1.200 RPM      | 800 RPM         |                 |                 |                 |                  |                 |                 |
|---------|--------------------|----------------|-----------------|-----------------|-----------------|-----------------|------------------|-----------------|-----------------|
|         | Ø<br>DIÂM.<br>(mm) | FL 2024<br>(A) | FHE 2324<br>(B) | FHE 3024<br>(C) | FHI 3324<br>(D) | FHI 3624<br>(E) | FHI 3824<br>(E+) | FHI 4024<br>(F) | FHI 4324<br>(G) |
| ILM38xx | 3,8                | •              | •               | •               | •               | •               | •                | •               | •               |
| ILM40xx | 4,0                | •              | •               | •               | •               | •               | •                | •               | •               |
| ILM45xx | 4,5                | •              | •               | •               | •               | •               | •                | •               | •               |

# ZYGOMATIC PLUS

## THE IMPLANTS

Zygomatic Plus implants are manufactured from Grade IV Cold Worked titanium and feature a dual acid-etching (DAA) surface treatment followed by the application of nano-hydroxyapatite (HAnano) in the threaded area, while the cervical area is treated only with HAnano without prior DAA.

Their insertion is performed using a carrier wrench for a contra-angle, which is provided in the specific surgical kit for this line. Installation is completed with a specific wrench for the Zygomatic implant, which can be attached to a hand wrench and/or torque wrench.


Zygomatic Plus implants must be placed in the upper maxillary arch to provide support for fixed dental prostheses in patients with fully edentulous maxillae. The Zygomatic Plus implant by S.I.N. has both threaded and non-threaded sections. The threaded portion of the implant, measuring 17.10 mm, is located in its apical part and is designed to provide anchorage in the patient's zygomatic bone. The surgeon should expect high-density bone in this area.

The smooth, non-threaded surface is intended to achieve passive stability over the alveolar process with a press-fit insertion into surgical sockets or resting against the bone, depending on the chosen technique. All implants are suitable for immediate loading when good primary stability is achieved and with appropriate occlusal loading.

The Zygomatic Plus line features a 16° Cone Morse connection (MT 16°) and comes in lengths of 30; 32.5; 35; 37.5; 40; 42.5; 45; 47.5; 50; 52.5; 55; 57.5; 60; and 62.5 mm, with a diameter of 4.0 mm, as detailed in the table below:

### TECHNICAL MEASUREMENTS

|            | A                         | B                       | C           |
|------------|---------------------------|-------------------------|-------------|
| CODE       | PLATFORM DIAMETER<br>(mm) | APICAL DIAMETER<br>(mm) | LENGTH (mm) |
| ILMZ 4030N | Ø4,0                      | Ø2,0                    | 30,0        |
| ILMZ 4032N | Ø4,0                      | Ø2,0                    | 32,5        |
| ILMZ 4035N | Ø4,0                      | Ø2,0                    | 35,0        |
| ILMZ 4037N | Ø4,0                      | Ø4,0                    | 37,5        |
| ILMZ 4040N | Ø4,0                      | Ø2,0                    | 40,0        |
| ILMZ 4042N | Ø4,0                      | Ø2,0                    | 42,5        |
| ILMZ 4045N | Ø4,0                      | Ø2,0                    | 45,0        |
| ILMZ 4047N | Ø4,0                      | Ø2,0                    | 47,5        |
| ILMZ 4050N | Ø4,0                      | Ø2,0                    | 50,0        |
| ILMZ 4052N | Ø4,0                      | Ø2,0                    | 52,5        |
| ILMZ 4055N | Ø4,0                      | Ø2,0                    | 55,0        |
| ILMZ 4057N | Ø4,0                      | Ø2,0                    | 57,5        |
| ILMZ 4060N | Ø4,0                      | Ø2,0                    | 60,0        |
| ILMZ 4062N | Ø4,0                      | Ø2,0                    | 62,5        |



## INSTALLATION


In severely resorbed maxillae, Zygomatic Plus implants are placed along with two or four conventional anterior implants. The Zygomatic Plus implants are anchored in the zygomatic bone. Typically, the prosthetic platform of the implant emerges slightly in the palatal position between the second premolar and the first upper molar region.

Two to four anterior implants should be installed axially in the palatal position, and two zygomatic implants should be tilted to provide an optimal prosthetic emergence, preferably at the center of the alveolus.

A 45° palatal incision along the entire bone crest, combined with a full-thickness flap from the maxillary crest to the zygomatic bone support and the identification of the infraorbital nerve, are the first steps of this surgery.

For the surgical procedure, if the patient has a complete lower denture, its removal may help manage the typically long instruments required for the placement of zygomatic implants.


Different drills with increasing diameters are used, concluding with the insertion of the ultra-threaded zygomatic implant at low speed. The implant length is selected using a probe and can range from 30.0 mm to 62.5 mm.



*Example of insertion of Zygomatic Plus implant associated with anterior implants.*



*Detail of the Zygomatic Plus implant and its anatomical relationship with the ocular cavity.*



*Infraorbital foramen.*


The technique for installation can be the one recommended by Branemark, with intrasinusal drilling, the Slot technique, where access and detachment of the maxillary sinus membrane are performed, or the exteriorized technique, where the implant will touch the wall of the maxillary sinus.

The most posterior implant is installed first; the palatine entry is made in the molar region, with the implant passing slightly behind the support and drilling the zygomatic bone from the medial face. The entry into the zygoma should be low and posterior, and the soft tissues should be carefully retracted. It is important to use sharp helical burs to avoid surgical accidents, as well as overheating.


The second implant is installed in the premolar region, following the lower zygomatic crest into the sinus and drilling the medial part of the zygomatic bone.

Ideally, each implant should be supported by the surrounding bone, both at the neck and at the apex. Often, the bone in the crest is very thin, so it is easy to create an entry hole larger than the implant. In patients with very severe resorption, the entry is sometimes located in the compact palatine bone. Therefore, it is important that the palatine hole and the zygoma hole have exactly the same direction to avoid stresses and difficulties during implant installation.

The maximum torque for installing the Zygomatic Plus, without causing damage or fractures, is 80 N/cm.



*Posterior margin of the zygomatic bone.*



*Anatomical references to protect the orbital cavity from injuries.*



*Viewing window.*

## KIT


For the installation of Zygomatic implants, it is necessary to use the specific kit for this line. A complete and compact kit with a linear and intuitive sequence.



| CODE      | DESCRIPTION                                |
|-----------|--------------------------------------------|
| FRLZ 27   | LAUNCH DRILL Ø2.35MM                       |
| FHZ 2030  | HELICOIDAL DRILL PILOTO Ø2.75MM            |
| FHZ 2932  | HELICOIDAL DRILL Ø3.0MM                    |
| FHZ 3234  | HELICOIDAL DRILL Ø3.2MM                    |
| FBD 40    | DIAMOND DRILL Ø4,0X40MM GROSSA             |
| FBD 40E   | SPHERICAL DIAMOND DRILL Ø4,0MM             |
| CQCA 27   | COUNTER ANGLE WRENCH 4.0MM                 |
| CTHA 1224 | COUNTER ANGLE WRENCH HEX.1.2MM EST. MEDIUM |
| FL 20M    | LAUNCH DRILL Ø2.0MM MEDIA                  |
| CCM 01L   | RATCHET KEY FOR IMPLANT CM LONG            |
| CCM 01M   | RATCHET KEY FOR IMPLANT CM MEDIUM          |
| CDH 1224  | DIGITAL WRENCH HEX 1.2 MEDIUM              |
| CMZ       | WRENCH FOR ZYGOMATIC IMPLANT               |
| SOPZ      | ZYGOMATIC DEPTH PROBE                      |

## MILLING

For implant installation, the following milling sequence should be followed.




| CODE      | REFERÊNCIA                             |
|-----------|----------------------------------------|
| FRLZ 27M  | LAUNCH DRILL Ø2.35MM MEDIUM            |
| FHZ 2030M | HELICOIDAL DRILL PILOTO Ø2.75MM MEDIUM |
| FHZ 2932M | HELICOIDAL DRILL Ø3.0MM MEDIUM         |
| FHZ 3234M | HELICOIDAL DRILL Ø3.2MM MEDIUM         |

\*The Launch and Helical Burs (A, B, C, and D) are also available in a short option to facilitate milling in specific cases. This version is sold separately.

## PROSTHETIC PLANNING

Angled abutments are recommended for tilted posterior implants. Install the mini-abutment with a torque of 20 N.cm using the Hex 1.2 wrench from the S.I.N. prosthetic kit.


The recommendation for axially installed implants is to use a straight mini-abutment with a torque of 20 N.cm using the wrench specifically for mini-abutments.



## CONVENTIONAL REHABILITATION (ANALOGUE)

Precisely place the impression transfer of the mini abutment onto the component and only tighten the screw manually or with the help of a 1.2 hexagonal wrench.

For taking the impression, it is recommended to use a standard elastomeric impression material (e.g., polyvinylsiloxane).



When the impression is ready, place the analogs and send the impression to the prosthetic laboratory. This final gypsum model will be used in the next steps of the rehabilitation process.



Wax the structure of the bar according to the available interocclusal space of the patient.

Place the bar and check its alignment with the model. Produce the final prosthesis based on the previously verified bar.

Install the definitive fixed prosthesis. Check that the final torque of the cylinder screws is 10 N.cm.

Check and adjust the correct occlusion.



# S.I.N. SOLUTIONS

The digital workflow is the sequence of steps necessary to automate processes, based on a well-defined set of rules that allow its transmission from one person to another.

Thanks to technological advancements, lower costs, and greater dissemination, more and more professionals are able to work with the digital workflow.

The digital workflow consists of an intraoral image, an impression, a model, or a digital wax-up to obtain a virtual model and enable planning and preparation of rehabilitation in computer software (CAD).

The software files are sent to a milling machine that produces or materializes the corresponding files (CAM).

STL stands for stereolithography. The format of an STL file uses a series of interconnected triangles to recreate the geometry of the solid model.

Due to the contrast achieved, colored scanning can help identify the tooth color and the finish line of the tooth preparation. Virtual models can also be used for complex prosthetic rehabilitations, orthodontic planning, osteotomies, surgical guides, and guided surgery.

• J Prosthet Dent. 2019 feb;121(2):246-251

***Accuracy of digital technologies for the scanning of facial, skeletal, and intraoral tissues: a systematic review.***

Bohner L, Gamba DD, Hanisch M, Marcio BS, Tortamano Neto P, Laganá DC, Sesma N.




## DIGITAL IMPRESSIONS

Specific transfer devices, called scanning JIGS, are used for intraoral scanning. These devices will later be converted into digital transfers, preventing the depth and amount of gingival tissue from interfering with the impression.

- Int J Oral Maxillofac Implants 2014;29(4):836-845

***Accuracy of implant impressions for partially and completely edentulous patients: a systematic review.***

Papaspyridakos P, Chen CJ, Gallucci GO, Doukoudakis A, Weber HP, Chronopoulos V.



Escaneo digital

## INTRAORAL SCANNER

Image capture is performed directly in the oral cavity. This eliminates clinical steps and avoids errors caused by distortions from impression materials and model fabrication.

This technique provides greater comfort for the patient, eliminating the need for a tray, impression materials, and gypsum models.

The intraoral scanner can be divided into two main units: the operational unit, which consists of the optical reader, and the processing control monitor, which guides the professional during image capture. Through the optical reader, light is incident on the patient's intraoral anatomical structures, which is reflected and detected by a sensor that is at a fixed distance and calculates the angle of reflection of the incident light beam through a process called triangulation.

• J Prosthet Dent. 2019 feb;121(2):246-251

***Accuracy of digital technologies for the scanning of facial, skeletal, and intraoral tissues: a systematic review.***

Bohner L, Gamba DD, Hanisch M, Marcio BS, Tortamano Neto P, Laganá DC, Sesma N.

• BMC Oral Health. 2020; 20:263-284

***Trueness of 12 intraoral scanners in the fullarch implant impression: a comparative in vitro study.***

Mangano FG., Admakin O., Bonacina M., Lerner H., Rutkunas V., Mangano C.



Intraoral scanner

## PROSTHETIC REHABILITATION

The provisional and definitive prostheses generated by the digital workflow can vary in terms of manufacturing types and materials.

Typically, they are milled using materials such as PMMA, zirconia, titanium, and others, or 3D printed using various types of resins specifically designed for use in dentistry.

- ***Precision and Accuracy of a Digital Impression Scanner in Full-Arch Implant Rehabilitation.***

Paolo Pesce, Francesco Pera, Paolo Setti, Maria Menini

Int J Prosthodont Mar/Apr 2018;31(2):171-175. doi: 10.11607/ijp.5535.

- ***Metal free, full arch, fixed prosthesis for edentulous mandible rehabilitation on four implants.***


Alfredo Passaretti 1, Giulia Petroni 1, Giovanna Miracolo 2, Valeria Savoia 2, Angelo Perpetuini 2, Andrea Cicconetti 3

J Prosthodont Res 2018 Apr;62(2):264-267. doi: 10.1016/j.jpor.2017.10.002. Epub 2017 Dec 7.


- ***Retrospective 2- to 7-Year Follow-Up Study of 20 Double Full-Arch Implant-Supported Monolithic Zirconia Fixed Prostheses: Measurements and Recommendations for Optimal Design.***

Fernando Rojas Vizcaya

Prosthodont 2018 Jul;27(6):501-508. doi: 10.1111/jopr.12528. Epub 2016 Aug 29.



Digital planning



Model with Jigs



Milled prosthesis



Milled prosthesis installed

# CONSIDERATIONS ON THE FULL ARCH CONCEPT

This concept is designed to make the most of the available bone and allow for immediate loading.

Points to consider when planning a full arch case:

- High primary stability
- Insertion torque: a minimum of 45 Ncm.
- If an insertion torque lower than this is recorded, we recommend a healing period before installing the provisional or definitive prosthesis.

If there are extraction sites, it is advisable to clean them before installing the implants.

It is advisable to place the implants between the post-extraction sockets.

To reduce cantilevers, tilt the posterior implants up to a maximum of 45°.

For tilted posterior implants, plan the screw access holes so that they are located in the occlusal plane of the first molar, second premolar, or first premolar.

# SCIENTIFIC PUBLICATIONS

- › Aires I, Berger J. Planning Implant Placement on 3D stereolithographic Models Applied with Immediate Loading of Implant-Supported Hybrid Protheses After Multiple Extractions: A Case Series. *Int J Oral Maxillofac Implants*. 2016; 31(1):172-178. doi: 10.11607/jomi.4186
- › Bohner L, Gamba DD, Hanisch M, Marcio BS, Tortamano Neto P, Laganá DC, Sesma N. Accuracy of digital technologies for the scanning of facial, skeletal, and intraoral tissues: a systematic review. *J Prosthet Dent*. 2019 feb;121(2):246-251
- › Bohner L, Gamba BB, Hanisch M, Marcio BS, Tortamano Neto P, Laganá DC, Sesma N. Comparison of conventional, photogrammetry, and intraoral scanning accuracy of complete arch implant impression procedures evaluated with a coordinate measuring machine. *Jprosthet Dent*. 2021 mar;125(3):470-478
- › Brånemark PI, Svensson B, van Steenberghe D. Ten-year survival rates of fixed protheses on four or six implants ad modum Brånemark in full edentulism. *Clin Oral Implants Res*. 1995 Dec;6(4):227-31
- › Carlos Aparicio, Carolina Manresa, Karen Francisco, Pedro Claros, Javier Al Andez, Oscar Gonzalez-Martin & Tomas Albrektsson. Zygomatic implants: indications, techniques and outcomes, and the Zygomatic Success Code. *Periodontology 2000*, Vol. 66, 2014, 41-58
- › Chrcanovic BR, Albrektsson T, Wennerberg A. Survival and Complications of Zygomatic Implants: An Updated Systematic Review. *J Oral Maxillofac Surg*. 2016 Oct;74(10):1949-64.
- › Duarte F, Ramos C, Silva JN. Immediate function with four zygomatic implants in patients with extreme maxillary atrophy – Case series. *Journal of Surgery, Periodontology and Implant Research* 2019; 1(2):51-55
- › Graves S, Mahler BA, Javid B, Armellini D, Jensen OT. Maxillary all-on-four therapy using angled implants: a 16-month study of 1110 implants in 276 jaws. *Dent Clin North Am* 2011;55:779-794.
- › Holtzclaw D. All-on-4® Implant Treatment: Common Pitfalls and Methods to Overcome Them. *Compend Contin Educ Dent*. 2016;37(7):458-465
- › Holtzclaw D, Telles R. Pterygoid Fixated Arch Stabilization Technique (PFAST): A Retrospective Study of Pterygoid Dental Implants used for Immediately Loaded Full Arch Prosthetics. *The Journal of Implant & Advanced Clinical Dentistry*. 2018 Dec;10(7):6-17
- › Horita S, Sugiura T, Yamamoto K, Murakami K, Imai Y, Kirita T. Biomechanical analysis of immediately loaded implants according to the “Allon-Four” concept. *J Prosthodont Res*. 2017 Apr;61(2):123-132.
- › Jensen OT, Adams MW. Secondary stabilization of maxillary m-4 treatment with unstable implants for immediate function: biomechanical considerations and report of 10 cases after 1 year in function. *Int J Oral Maxillofac Implants*. 2014;29(2):232-40.

- › Lombardo G, D'Agostino A, Trevisiol L, Romanelli MG, Mascellaro A, Gomez-Lira M, Pardo A, Favero V, Nocini PF. Clinical, microbiologic and radiologic assessment of soft and hard tissues surroundingzygomatic implants: a retrospective study. *Oral Surg Oral Med Oral Pathol Oral Radiol*. 2016 Nov;122(5):537-546
- › Ma B., Yue X., Sun Y., Peng L., Geng W. Accuracy of photogrammetry, intraoral scanning, and conventional impression techniques for complete-arch implant rehabilitation: an in vitro comparative study. *Bmc Oral Health*. 2021; 21:636-644
- › Malchiodi L, Moro T, Cattina DP, Cucchi A, Ghensi P, Nocini PF. Implant rehabilitation of the edentulous jaws: Does tilting of posterior implants at an angle greater than 45 affect bone resorption and implant success?: A retrospective study. *Clin Implant Dent Relat Res* 2018 Oct;20(5):867-874
- › Mangano FG., Admakin O., Bonacina M., Lerner H., Rutkunas V., Mangano C. Trueness of 12 intraoral scanners in the fullarch implant impression: a comparative in vitro study. *Bmc Oral Health*. 2020; 20:263-284
- › Marcelo de Carvalho, Liliane Pacheco de Carvalho, Rogerio de Lima Romeiro, Carlos Eduardo Francischone, Bruno Salles Sotto-Maior, Fabio Bezerra. Nova Proposta Para Reabilitação de Maxila Atrófica: Implante Inclinado LONG. *Int J Oral Maxillofac Implants* 2016;31:1017-1022. (Portuguese Edition)
- › Mattsson T, Köndell PA, Gynther GW, Fredholm U, Bolin A. Implant treatment without bone grafting in severely resorbed edentulous maxillae. *J Oral Maxillofac Surg* 1999; 57(3):281-287
- › Maló P, de Araújo Nobre M, Lopes A, Moss SM, Molina GJ. A longitudinal study of the survival of Allon- 4 implants in the mandible with up to 10 years of follow-up. *J Am Dent Assoc*. 2011;142(3):310-20.
- › Maló P, de Araújo Nobre M, Lopes A, Ferro A, Gravito I. All-on-4® Treatment Concept for the Rehabilitation of the Completely Edentulous Mandible: A 7-Year Clinical and 5-Year Radiographic Retrospective Case Series with Risk Assessment for Implant Failure and Marginal Bone Level. *Clin Implant Dent Relat Res*. 2015;17(2):531-41.
- › Maló P, Rangert B., Nobre M. All-on-4 immediate-function concept with Brånemark System implants for completely edentulous maxillae: a 1-year retrospective clinical study. *Clin Implant Dent Relat Res*. 2005;7 Suppl 1:S88-94
- › Martelli N, Serrano C, Van den Brink H, Pineau J, Prognon P, Borget.I. Advantages and disadvantages of 3-dimensional printing in surgery: A systematic review. *Surgery*. 2016 June; 159(6):1485-1500. doi: 10.1016/j.surg.2015.12.017.
- › Meloni SM, De Riu G, Pisano M, Lolli FM, Deledda A, Guglielmo Campus G, et al. Implant Restoration of Edentulous Jaws with 3D Software Planning, Guided Surgery, Immediate Loading, and CAD-CAM Full Arch Frameworks. *Int J Dent*. 2013;2013:683423. doi: 10.1155/2013/683423.
- › Padhye OV, Herekar M, Patil V, Mulani S, Sethi M, Fernandes A. Stress Distribution in Bone and Implants in Mandibular 6-Implant Supported Cantilevered Fixed Prosthesis: A 3D Finite Element Study *Implant Dent*. 2015 Dec;24(6):680-5.

- › Patzelt SB, Bahat O, Reynolds MA, Strub JR. The all-on-four treatment concept: a systematic review. *Clin Implant Dent Relat Res.* 2014 Dec;16(6):836-55
- › Peñarrocha M, Carrillo C, Boronat A, Balaguer J, Peñarrocha M. Palatal Positioning of Implants in Severely Resorbed Edentulous Maxillae. *Int J Oral Maxillofac Implants* 2009 May-Jun;24(3):527-33
- › Papaspyridakos P, Chen CJ, Gallucci GO, Doukoudakis A, Weber HP, Chronopoulos V. Accuracy of implant impressions for partially and completely edentulous patients: a systematic review. *Int J Oral Maxillofac Implants* 2014;29(4):836-845
- › Ozan O, Turkyilmaz I, Ersoy AE, McGlumphy EA, Rosenstiel SF. Clinical accuracy of 3 different types of computed tomography-derived stereolithographic surgical guides in implant placement. *J Oral Maxillofac Surg.* 2009 Feb;67(2):394-401. doi: 10.1016/j.joms.2008.09.033.
- › Romanos GE, Nentwig GH. Immediate functional loading in the maxilla using implants with platform switching: five-year results. *Int J Oral Maxillofac Implants* 2009;24:1106–1112.
- › Reiz SD, Neugebauer J, Karapetian VE, Ritter L. Cerec Meets Galileos - Integrated Implantology for Completely Virtual Implant Planning. *Int J Comput Dent.* 2014 17(2):145-157.
- › Ujigawa K, Kato Y, Kizu Y, Tonogi M, Yamane GY. Three – dimensional finite elemental analysis of zygomatic implants in craniofacial structures. *Int J Oral Maxillofac Surg.* 2007 Jul;36(7):620-5

# OUR GLOBAL PRESENCE



POINT YOUR PHONE'S CAMERA AT THE  
QR CODE AND SEE WHERE S.I.N. IS PRESENT

## HEADQUARTERS

2140 Vereador Abel Ferreira Av  
Jardim Anália Franco  
São Paulo – SP - Brazil

## FACTORY

421 Soldado Ocimar Guimarães  
da Silva St - Jardim Anália Franco  
São Paulo – SP - Brazil

## PORTUGAL BRANCH

General Ferreira Martins St, 10 8D -  
1495-137 Algés - Portugal  
+351 21 412-0336

[sinimplantsystem.com](http://sinimplantsystem.com)





Get to know Implantat, the educational streaming of S.I.N.

 [implantat.global](https://implantat.global)

[internacional@sinimplantsystem.com](mailto:internacional@sinimplantsystem.com)  
[www.sinimplantsystem.com](https://www.sinimplantsystem.com)

---

Visit our social networks



[@sinimplantglobal](https://www.instagram.com/sinimplantglobal)



[@sinimplantglobal](https://www.facebook.com/sinimplantglobal)



[S.I.N. Implant System](https://www.youtube.com/S.I.N.ImplantSystem)



[/sin\\_implant](https://www.tiktok.com/@sin_implant)