Dr. Fábio J. B. Bezerra

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- Sohrabi K, Mushantat A, Esfandiari S, Feine J. How successful are small-diameter implants? A literature review. *Clin Oral Implants Res.* 2012; 23 (5): 515-25.
- Gaviria L, Salcido JP, Guda T, Ong JL. Current trends in dental implants. *J Korean Assoc Oral Maxillofac Surg.* 2014; 40 (2): 50-60.
- Martinez EF, Ishikawa GJ, de Lemos AB, Barbosa Bezerra FJ, Speardio M, Napimoga MH. Evaluation of a Titanium Surface Treated with Hydroxyapatite Nanocrystals on Osteoblastic Cell Behavior: An In Vitro Study. *Int J Oral Maxillofac Implants.* 2018; 33 (3): 597-602.
- Rossi MC, Bezerra FJB, Silva RA, Crulhas BP, Fernandes CJC, Nascimento AS, et al. Titanium-released from dental implant enhances pre-osteoblast adhesion by ROS modulating crucial intracellular pathways. *J Biomed Mater Res A.* 2017; 105 (11): 2968-76.
- Fernandes CJC, Bezerra F, do Carmo MDD, Feltran GS, Rossi MC, da Silva RA, et al. CoCr-enriched medium modulates integrin-based downstream signaling and requires a set of inflammatory genes reprogramming in vitro. *J Biomed Mater Res A.* 2018; 106 (3): 839-49.
- Di Liddo R, Paganin P, Lora S, Dalzoppo D, Giraudo C, Miotto D, et al. Poly-epsilon-caprolactone composite scaffolds for bone repair. *Int J Mol Med.* 2014; 34 (6): 1537-46.
- Salou L, Hoornaert A, Stanovici J, Briand S, Louarn G, Layrolle P. Comparative bone tissue integration of nanostructured and microroughened dental implants. *Nanomedicine (Lond).* 2015; 10 (5): 741-51.
- Wennerberg A, Albrektsson T. Effects of titanium surface topography on bone integration: a systematic review. *Clin Oral Implants Res.* 2009; 20 Suppl 4: 172-84.
- Rosa MB, Albrektsson T, Francischone CE, Filho HO, Wennerberg A. Micrometric characterization of the implant surfaces from the five largest companies in Brazil, the second largest worldwide implant market. *Int J Oral Maxillofac Implants.* 2013; 28 (2): 358-65.
- Zambuzzi WF, Bonfante EA, Jimbo R, Hayashi M, Andersson M, Alves G, et al. Nanometer scale titanium surface texturing are detected by signaling pathways involving transient FAK and Src activations. *PloS one.* 2014; 9 (7): e95662.
- Da Silva RA, da Silva Feltran G, Ferreira MR, Wood PF, Bezerra F, Zambuzzi WF. The Impact of Bioactive Surfaces in the Early Stages of Osseointegration: An In Vitro Comparative Study Evaluating the HAnano(R) and SLActive(R) Super Hydrophilic Surfaces. *Biomed Res Int.* 2020; 2020: 3026893.
- Bezerra F, Ferreira MR, Fontes GN, da Costa Fernandes CJ, Andia DC, Cruz NC, et al. Nano hydroxyapatite-blasted titanium surface affects pre-osteoblast morphology by modulating critical intracellular pathways. *Biotechnol Bioeng.* 2017; 114 (8): 1888-98.
- Barkarmo S, Andersson M, Currie F, Kjellin P, Jimbo R, Johansson CB, et al. Enhanced bone healing around nanohydroxyapatite-coated polyetheretherketone implants: An experimental study in rabbit bone. *Journal of biomaterials applications.* 2014; 29 (5): 737-47.
- Jimbo R, Sotres J, Johansson C, Breiding K, Currie F, Wennerberg A. The biological response to three different nanostructures applied on smooth implant surfaces. *Clin Oral Implants Res.* 2012; 23 (6): 706-12.
- Agarwal R, Garcia AJ. Biomaterial strategies for engineering implants for enhanced osseointegration and bone repair. *Adv Drug Deliv Rev.* 2015; 94: 53-62.
- Gemini-Piperni S, Milani R, Bertazzo S, Peppelenbosch M, Takamori ER, Granjeiro JM, et al. Kinome profiling of osteoblasts on hydroxyapatite opens new avenues on biomaterial cell signaling. *Biotechnol Bioeng.* 2014; 111 (9): 1900-5.
- Gemini-Piperni S, Takamori ER, Sartoretto SC, Paiva KB, Granjeiro JM, de Oliveira RC, et al. Cellular behavior as a dynamic field for exploring bone bioengineering: a closer look at cell-biomaterial interface. *Arch Biochem Biophys.* 2014; 561: 88-98.
- Henkel J, Woodruff MA, Epari DR, Steck R, Glatt V, Dickinson IC, et al. Bone Regeneration Based on Tissue Engineering Conceptions - A 21st Century Perspective. *Bone Res.* 2013; 1 (3): 216-48.
- Bertazzo S, Zambuzzi WF, Campos DD, Ferreira CV, Bertran CA. A simple method for enhancing cell adhesion to hydroxyapatite surface. *Clin Oral Implants Res.* 2010; 21 (12): 1411-3.
- Bertazzo S, Zambuzzi WF, Campos DD, Ogeda TL, Ferreira CV, Bertran CA. Hydroxyapatite surface solubility and effect on cell adhesion. *Colloids Surf B Biointerfaces.* 2010; 78 (2): 177-84.
- Zambuzzi WF, Fernandes GV, Iano FG, Fernandes Mda S, Granjeiro JM, Oliveira RC. Exploring anorganic bovine bone granules as osteoblast carriers for bone bioengineering: a study in rat critical-size calvarial defects. *Brazilian dental journal.* 2012; 23 (4): 315-21.
- Zambuzzi WF, Ferreira CV, Granjeiro JM, Aoyama H. Biological behavior of pre-osteoblasts on natural hydroxyapatite: a study of signaling molecules from attachment to differentiation. *J Biomed Mater Res A.* 2011; 97 (2): 193-200.
- Zambuzzi WF, Coelho PG, Alves GG, Granjeiro JM. Intracellular signal transduction as a factor in the development of "smart" biomaterials for bone tissue engineering. *Biotechnol Bioeng.* 2011; 108 (6): 1246-50.
- Mello-Machado RC, Sartoretto SC, Granjeiro JM, Calasans-Maia JA, de Uzeda M, Mourao C, et al. Osseodensification enables bone healing chambers with improved low-density bone site primary stability: an in vivo study. *Sci Rep.* 2021; 11 (1): 15436.
- Sartoretto SC, Calasans-Maia J, Resende R, Camara E, Ghiraldini B, Barbosa Bezerra FJ, et al. The Influence of Nanostructured Hydroxyapatite Surface in the Early Stages of Osseointegration: A Multiparameter Animal Study in Low-Density Bone. *Int J Nanomedicine.* 2020; 15: 8803-17.
- Hu N, Feng C, Jiang Y, Miao Q, Liu H. Regulative Effect of Mir-205 on Osteogenic Differentiation of Bone Mesenchymal Stem Cells (BMSCs): Possible Role of SATB2/Runx2 and ERK/MAPK Pathway. *Int J Mol Sci.* 2015; 16 (5): 10491-506.
- Jiang T, Guo L, Ni S, Zhao Y. Upregulation of cell proliferation via Shc and ERK1/2 MAPK signaling in SaOS-2 osteoblasts grown on magnesium alloy surface coating with tricalcium phosphate. *J Mater Sci Mater Med.* 2015; 26 (4): 158.
- Asati V, Mahapatra DK, Bharti SK. PI3K/Akt/mTOR and Ras/Raf/MEK/ERK signaling pathways inhibitors as anticancer agents: Structural and pharmacological perspectives. *Eur J Med Chem.* 2016; 109: 314-41.
- Branemark PI, Adell R, Breine U, Hansson BO, Lindstrom J, Ohlsson A. Intra-osseous anchorage of dental prostheses. I. Experimental studies. *Scand J Plast Reconstr Surg.* 1969; 3 (2): 81-100.
- Albrektsson T, Branemark PI, Hansson HA, Lindstrom J. Osseointegrated titanium implants. Requirements for ensuring a long-lasting, direct bone-to-implant anchorage in man. *Acta Orthop Scand.* 1981; 52 (2): 155-70.
- Wennerberg A, Albrektsson T, Johansson C, Andersson B. Experimental study of turned and grit-blasted screw-shaped implants with special emphasis on effects of blasting material and surface topography. *Biomaterials.* 1996; 17 (1): 15-22.
- Coelho PG, Lemons JE. Physico/chemical characterization and in vivo evaluation of nanothickness bioceramic depositions on alumina-blasted/acid-etched Ti-6Al-4V implant surfaces. *J Biomed Mater Res A.* 2009; 90 (2): 351-61.
- Coelho PG, Jimbo R, Tovar N, Bonfante EA. Osseointegration: hierarchical designing encompassing the micrometer, micrometer, and nanometer length scales. *Dent Mater.* 2015; 31 (1): 37-52.
- Albrektsson T, Wennerberg A. Oral implant surfaces: Part 2--review focusing on clinical knowledge of different surfaces. *Int J Prosthodont.* 2004; 17 (5): 544-64.
- Yoon SR, Cha JK, Lim HC, Lee JS, Choi SH, Jung UW. De novo bone formation underneath the sinus membrane supported by a bo-

- ne patch: a pilot experiment in rabbit sinus model. *Clin Oral Implants Res.* 2016.
36. **Buser D, Brogini N, Wieland M, Schenk RK, Denzer AJ, Cochran DL, et al.** Enhanced bone apposition to a chemically modified SLA titanium surface. *J Dent Res.* 2004; 83 (7): 529-33.
 37. **Cochran DL, Schenk RK, Lussi A, Higginbottom FL, Buser D.** Bone response to unloaded and loaded titanium implants with a sandblasted and acid-etched surface: a histometric study in the canine mandible. *J Biomed Mater Res.* 1998; 40 (1): 1-11.
 38. **Pan J, Liao H, Leygraf C, Thierry D, Li J.** Variation of oxide films on titanium induced by osteoblast-like cell culture and the influence of an H₂O₂ pretreatment. *J Biomed Mater Res.* 1998; 40 (2): 244-56.
 39. **Nanci A, Wuest JD, Peru L, Brunet P, Sharma V, Zalzal S, et al.** Chemical modification of titanium surfaces for covalent attachment of biological molecules. *J Biomed Mater Res.* 1998; 40 (2): 324-35.
 40. **De Oliveira PT, Nanci A.** Nanotexturing of titanium-based surfaces upregulates expression of bone sialoprotein and osteopontin by cultured osteogenic cells. *Biomaterials.* 2004; 25 (3): 403-13.
 41. **Shi GS, Ren LF, Wang LZ, Lin HS, Wang SB, Tong YQ.** H₂O₂/HCl and heat-treated Ti-6Al-4V stimulates pre-osteoblast proliferation and differentiation. *Oral Surg Oral Med Oral Pathol Oral Radiol Endod.* 2009; 108 (3): 368-75.
 42. **Dohan Ehrenfest DM, Coelho PG, Kang BS, Sul YT, Albrektsson T.** Classification of osseointegrated implant surfaces: materials, chemistry and topography. *Trends Biotechnol.* 2010; 28 (4): 198-206.
 43. **Oh SH, Finones RR, Daraio C, Chen LH, Jin S.** Growth of nanoscale hydroxyapatite using chemically treated titanium oxide nanotubes. *Biomaterials.* 2005; 26 (24): 4938-43.
 44. **Hall J, Miranda-Burgos P, Sennerby L.** Stimulation of directed bone growth at oxidized titanium implants by macroscopic grooves: an in vivo study. *Clin Implant Dent Relat Res.* 2005; 7 Suppl 1: S76-82.
 45. **Weiner S, Traub W.** Bone structure: from angstroms to microns. *FASEB J.* 1992; 6 (3): 879-85.
 46. **Aumailley M, Gayraud B.** Structure and biological activity of the extracellular matrix. *J Mol Med (Berl).* 1998; 76 (3-4): 253-65.
 47. **Ellingsen JE, Johansson CB, Wennerberg A, Holmen A.** Improved retention and bone-to-implant contact with fluoride-modified titanium implants. *Int J Oral Maxillofac Implants.* 2004; 19 (5): 659-66.
 48. **Coelho PG, Freire JN, Granato R, Marin C, Bonfante EA, Gil JN, et al.** Bone mineral apposition rates at early implantation times around differently prepared titanium surfaces: a study in beagle dogs. *Int J Oral Maxillofac Implants.* 2011; 26 (1): 63-9.
 49. **Webster TJ, Ahn ES.** Nanostructured biomaterials for tissue engineering bone. *Adv Biochem Eng Biotechnol.* 2007; 103: 275-308.
 50. **Meirelles L, Arvidsson A, Andersson M, Kjellin P, Albrektsson T, Wennerberg A.** Nano hydroxyapatite structures influence early bone formation. *J Biomed Mater Res A.* 2008; 87 (2): 299-307.
 51. **Jimbo R, Xue Y, Hayashi M, Schwartz-Filho HO, Andersson M, Mustafa K, et al.** Genetic responses to nanostructured calcium-phosphate-coated implants. *J Dent Res.* 2011; 90 (12): 1422-7.
 52. **Shunmugasamy VC, Gupta N, Pessoa RS, Janal MN, Coelho PG.** Influence of clinically relevant factors on the immediate biomechanical surrounding for a series of dental implant designs. *J Biomech Eng.* 2011; 133 (3): 031005.
 53. **Dang Y, Zhang L, Song W, Chang B, Han T, Zhang Y, et al.** In vivo osseointegration of Ti implants with a strontium-containing nanotubular coating. *Int J Nanomedicine.* 2016; 11: 1003-11.
 54. **Kulkarni M, Mazare A, Gongadze E, Perutkova S, Kralj-Iglic V, Milosev I, et al.** Titanium nanostructures for biomedical applications. *Nanotechnology.* 2015; 26 (6): 062002.
 55. **Longo G, Ioannidu CA, Scotto d'Abusco A, Superti F, Misiano C, Zanoni R, et al.** Improving Osteoblast Response In Vitro by a Nanostructured Thin Film with Titanium Carbide and Titanium Oxides Clustered around Graphitic Carbon. *PloS one.* 2016; 11 (3): e0152566.
 56. **Zambuzzi WF, Milani R, Teti A.** Expanding the role of Src and protein-tyrosine phosphatases balance in modulating osteoblast metabolism: lessons from mice. *Biochimie.* 2010; 92 (4): 327-32.
 57. **Zambuzzi WF, Bruni-Cardoso A, Granjeiro JM, Peppelenbosch MP, de Carvalho HF, Aoyama H, et al.** On the road to understanding of the osteoblast adhesion: cytoskeleton organization is rearranged by distinct signaling pathways. *J Cell Biochem.* 2009; 108 (1): 134-44.
 58. **Gottlander M, Johansson CB, Wennerberg A, Albrektsson T, Radin S, Ducheyne P.** Bone tissue reactions to an electrophoretically applied calcium phosphate coating. *Biomaterials.* 1997; 18 (7): 551-7.
 59. **Jimbo R, Coelho PG, Bryington M, Baldassarri M, Tovar N, Currie F, et al.** Nano hydroxyapatite-coated implants improve bone nanomechanical properties. *J Dent Res.* 2012; 91 (12): 1172-7.
 60. **Oliveira P, Coelho PG, Bergamo ETP, Witek L, Borges CA, Bezerra FB, et al.** Histological and Nanomechanical Properties of a New Nanometric Hydroxyapatite Implant Surface. An In Vivo Study in Diabetic Rats. *Materials (Basel).* 2020;13(24).
 61. **Oliveira P, de Melo Soares MS, Silveira ESAMM, Taba M, Jr., Pallioto DB, Messori MR, et al.** Influence of nano-hydroxyapatite coating implants on gene expression of osteogenic markers and micro-CT parameters. An in vivo study in diabetic rats. *J Biomed Mater Res A.* 2021; 109 (5): 682-94.