



Epikut


The new definition of epic

 S.I.N.



Epikut^S THE NEW DEFINITION OF EPIC

We have redefined the concept of the epic with EPIKUT. With a sharp design, compressive structure, and inverted double-support threads, this line offers greater clinical practicality, predictability, and high primary stability for those seeking superior results.



› Exclusive HAnano Plus Surface

Developed at leading universities in Sweden, this nanosurface accelerates osseointegration and promotes superior bone quality, backed by more than 50 preclinical studies.



› Apex

Stability and support for cases with low bone density.



› Exclusive Cervical Microthreads

Greater bone contact area and improved dissipation of occlusal forces.




› One Implant, Multiple Possibilities

Morse Taper Connection 11.5° and External Hexagon for easier clinical workflows.


› Double Acid Etching

Across the entire surface for Morse Taper and up to the cervical region in HE implants.



› Inverted Double-Support Threads

Ensure greater primary stability and insertion torque.




Epikut^S THE POWER OF BEING EPIC

Developed to enhance clinical practice, the Epikut S combines the unique features of the Epikut line with the 16° Morse Taper connection. The Epikut S offers predictability, high performance, and superior results. Features for those who want to be epic every day!




› Inverted Double-Support Threads

Ensure greater primary stability and insertion torque.



› More Options in the Portfolio

16° Morse Taper Connection to simplify clinical workflows.




› Exclusive Cervical Microthreads

Larger bone contact area and improved dissipation of occlusal forces.



› Apex


Stability and support for cases with low bone density.




› Accelerated Bone Healing

The high hydrophilicity, generated by an ultrathin and homogeneous layer of hydroxyapatite, enhances the activity of proteins involved in the osseointegration process.


★ EXTERNAL HEXAGONO




○ MORSE TAPER 11,5°



○ MORSE TAPER 16°





Epikut LONG


AN **EPIC** IMPLANT FOR
CASES WHERE YOU NEED
TO GO **BEYOND**

With a sharp, compressive design and inverted double-support threads, the Epikut Long line offers more size options to simplify your clinical workflow.



► An Implant with Many Possibilities

The Epikut Long line offers four different lengths (18, 20, 22, and 24 mm), ideal for full maxillary rehabilitation cases.




► Agility and Predictability

It allows for the rehabilitation of edentulous patients in a simpler, faster, and more predictable way, offering the possibility to avoid grafts or sinus lifts.

► Fit Precision

With unique abutments highly resistant to tension.




- Alta previsibilidad para la carga inmediata.
- Excelente anclaje en maxilas atroficas.
- Resolución del caso en un único procedimiento quirúrgico.




► Options for your preference:

Long implants are available in the Epikut line with a 11.5° Morse Taper connection and in the Epikut S line with a 16° Morse Taper connection.

MORSE TAPER LONG 16°




MORSE TAPER LONG 11,5°



EPIKUT SURGICAL KIT


MAXIMUM FUNCTIONALITY AND SIMPLICITY FOR YOUR SURGERIES



PRODUCT CODE: KCSE 02
ORGANIZING BOX CODE: COSE 02

EPIKUT LONG SURGICAL KIT

MAXIMUM FUNCTIONALITY AND SIMPLICITY FOR YOUR SURGERIES



CODE: KCSEL 02
ORGANIZING BOX CODE: COSEL 02

EPIKUT SAFE DRILL KIT

MAKING YOUR SURGERIES MORE PRACTICAL AND PRECISE

*The Epikut Safe Drill Kit is not compatible with the Epikut Long Surgical Kit.



For the installation of the Morse Taper to occur as recommended (subosseous), it is necessary to use a limiter 1.5 mm greater than the desired depth.

CODE: KESD 02
ORGANIZING BOX CODE: COESD 02

*Check product availability in your country.

*Check product availability in your country.



Get to know Implantat, the educational streaming of S.I.N.

 [implantat.global](https://www.implantat.global)

internacional@sinimplantsystem.com
www.sinimplantsystem.com

Visit our social networks



[@sinimplantglobal](https://www.instagram.com/sinimplantglobal)



[@sinimplantglobal](https://www.facebook.com/sinimplantglobal)



[S.I.N. Implant System](https://www.youtube.com/S.I.N.ImplantSystem)



[/sin_implant](https://www.tiktok.com/sin_implant)